



propper
Manufacturing Co., Inc.

36-04 Skillman Avenue
Long Island City, NY 11101
Tel (800) 832- 4300
Fax (718) 482-8909

www.propermfg.com

OK[®] Sterilizer Bags

Reorder #: 02100200, 02100300, 02101100

Quick Overview

These all-paper, high quality, gusseted sterilization packaging bags with pre-printed indicators are available in three sizes and are an ideal solution in rigid containers.

Details

Sterilization Paper Bag Construction

The exterior of each bag carries a pre-printed process indicator to allow easy identification of packs that have been subjected to the steam sterilization process. Superior quality medical grade paper and seals maintain the bags integrity even in the case of wet steam. The vertical seam is double bonded for security and there are designated areas to log the date, the operators' initials, and the contents.

Autoclave Tape Sealing

The bags are designed to be closed by folding over the top edge a number of times and then sealing this with autoclave tape. Their "fold and tape" design is ideal for easy and quick use. It also allows the prevention of the problems associated with imperfect heat sealing that may be due to faulty equipment or user error.

Manufacturer Recommendations for Use

Propper recommends OK Sterilizer Bags for use within wrapped sets or containment devices as a better alternative to paper-plastic pouches. All-paper bags allow for adequate air removal, steam contact, and drying due to steam penetration through all sides of the bag. For that purpose, paper-plastic pouches, according to the section "Paper plastic pouches" of ANSI/AAMI/ ST-79:2013, are not appropriate.

Features:

- Economical: The most economical packaging material for sterilization
- Surfaces have extra-smooth, machine-glazed paper finish to eliminate line
- All-paper bag will not become dry or brittle allowing each item to be stored for the maximum recommended period of time
- Each bag has a printed chemical indicator that turns brown after a successful steam sterilization process
- All sizes possess a gusset that allows bulky objects to be placed inside the bags.

*Each bag may be processed only once.