

## Instructions for Use: OK-Sonic™

### Product Description:

OK-Sonic™ is designed to monitor typical cleaning cycles of ultrasonic baths used in clinical and laboratory settings. The OK-Sonic™ monitor can be used for daily monitoring of cleaning efficacy of ultrasonic cleaner for unit release, or periodically for load release. OK-Sonic™ features two soil spots which will withstand water exposure without cavitation, and has been validated against traditional foil tests.

To use OK-Sonic™, slide an unprocessed indicator in the Monitor Strip Holder and hang monitor on rack within ultrasonic bath. Start your ultrasonic bath cycle as normal. At the completion of the cycle, remove the indicator from the holder and inspect for the presence of remaining blue test soil. Avoid contact with any remaining soil when the indicator is still wet. If there is any remaining soil, your test has failed; which should result in reviewing the cleaning cycle process and possible maintenance of your ultrasonic bath. If the indicator is free of any of remaining residue, your test has passed and it is safe to move forward to final confirmation of cleanliness.

### Instructions For Use:

1. Fill ultrasonic cleaner with appropriate cleaning solution. (Recommended concentration and temperature of cleaning solution depends on detergent; refer to detergent manufacturer's IFU for concentration and temperature of cleaning solution).
2. Degas the solution until solution is free of bubbles. (Recommended degassing time depends on ultrasonic unit; refer to ultrasonic manufacturer's IFU for degassing instructions).
3. Place OK-Sonic™ indicator into holder.
4. Place holder into ultrasonic cleaner in area to be tested.
5. Run cleaning cycle in accordance with ultrasonic manufacturer's instructions. At end of cycle, remove indicator from holder and check for presence of blue soil on indicator strip. Remaining blue soil indicates a fail.

**Shelf Life:** 24 months

**Storage prior to use:** Store at max temperature 30°C (86°F) and 30-70% RH

**Storage after use:** Store at 35°C (95°F) and 35-70% RH, OK-Sonic™ can be kept for long term record keeping

Reorder Number: 26210200

