

Propper Manufacturing Company

36-04 Skillman Avenue Long Island City, NY 11101

Tel. 718-392-6650

36-04 Skillman Avenue Long Island City, NY 11101 Tel. (718) 392-6650 • Fax. (718) 482-8909 www.proppermfg.com

March 25, 2024
Dear Customer,
Duo-Spore® Biological Indicator Strips (reorder # 26906500, 26906700, 26909200, 26909300, 26909400, 26909500, 26909600, 26909700, 26910600) have been deemed non-hazardous under the guidelines provided by OSHA Hazard Communication Standard 29 CFR 1910.1200; because of this, no Safety Data Sheet (SDS) is required for these products.
These products contain no known ecological hazards. The microorganisms in this product are non-pathogenic and are classified as Group 1 organisms by the World Health Organization (WHO). Products do not contain flammable liquids or gases. Standard transport conditions apply; no special conditions are required. Incineration or autoclaving at 121°C for 30 minutes are the recommended methods of disposal.
The adhesive present in the glassine packaging contains natural rubber latex.
Please contact Propper Manufacturing Company with any questions.
Sincerely,