

Infectious Waste Control

Biohazard & Decontamination Bags

Confidence

Decontamination and biohazard bags have pre-printed steam sterilization indicators that change from tan to black to indicate that the bag has been autoclaved

Visibility

Both bags display a highly visible biohazard symbol which are easily identified by hospital staff

Strength

Welded seams and reinforced corners provide durability assurance and a low risk of leakage

Biohazard Bags

Propper Orange Biohazard Bags are suitable for the collection, sterilization, and disposal of contaminated waste.



- Biohazard bags are made of high strength, 38µm thick polypropylene
- Propper biohazard bags are capable of withstanding high autoclave temperatures up to 286°F (141°C)
- Twist ties are supplied for easy closure of bags

Reorder No.	Qty/Case	Inches	Centimeters
32603500	200	14x19	35.6x48.3
32602300	200	19x23	8.3x58.4
32601900	200	25x35	63.5x88.9

Decontamination Bags

Propper Transparent Decontamination Bags are designed for the collection, transportation, and decontamination of medical or surgical equipment or waste. Transparency allows for easy identification of the contents of the bag. The special plastic construction enables immediate autoclaving when necessary.

- Propper Decontamination Bags are made from 51µm thick polypropylene
- Bags are capable of withstanding autoclave temperatures up to 286°F (141°C)
- Twist ties are supplied for easy closure of bags

Reorder No.	Qty/Case	Inches	Centimeters
32202400	200	25x35	63.5x88.9
32203000	200	24x30	48.26x58.42
32203600	200	14x19	35.56x48.26

