



Contact us to order Steam Chex today!*

customerservice@proppermfg.com

Tel. (718) 392-6650

Reorder No.: Quantity:

26810800 48 Rolls • $\frac{3}{4}$ " x 60 YDS

26810900 36 Rolls • 1" x 60 YDS


Blue Autoclave Indicator Tape and Steam Chex are pressure-sensitive process indicator tapes designed to seal packs that are exposed to steam sterilization.

Blue Autoclave Indicator Tape and Steam Chex are made of a saturated crepe paper printed with non-lead based ink, and coated with solvent free adhesive.

Blue Autoclave Indicator Tape and Steam Chex are not made with lead based ink or natural latex rubber.

Blue Autoclave Indicator Tape features a **WHITE** to dark brown/black color change

Steam Chex features a **GREEN** to dark brown/black color change



*Blue Autoclave Indicator Tape will no longer be available for purchase in the Fall of 2021.