



FDA Approved

Digital Thermometer

Oral or Armpit Measurement

Fast and Accurate Results

Easy to Clean and Disinfect

Easy-to-Read Digital Display

**Temperature Displays
in less than 1 minute**



Features

Reorder Number:
26915500

- ◆ **Ultrathin Thermometer**
- ◆ **Use 2 LR/sr-41 1.5V batteries**
- ◆ **Auto power-off after 10 minutes**

Product Specifications

Measurement Units : Celsius (°C)/Fahrenheit (°F)

Measurement Range : 90°F - 109.9°F

Measurement Accuracy : $\pm 0.2^{\circ}\text{F}$

Dimensions : 130mm x 19mm x 11mm

Weight : 10g