



April 11, 2016

Dear Customer,

Duo-Spore Biological Indicator Strips (reorder # 269065, 269067, 269092, 269093, 269094, 269095, 269096, 269097, 269106) have been deemed non-hazardous under the guidelines provided by OSHA Hazard Communication Standard 29 CFR 1910.1200; because of this, no Safety Data Sheet (SDS) is required for these products.

These products contain no known ecological hazards. The microorganisms in this product are non-pathogenic and are classified as Group I organisms by the World Health Organization (WHO). Products do not contain flammable liquids or gases. Standard transport conditions apply; no special conditions are required. Incineration or autoclaving at 121°C for 30 minutes are the recommended methods of disposal.

The adhesive present in the blue glassine packaging contains natural rubber latex.

Please contact Propper Manufacturing Company with any questions.

Sincerely,

Propper Manufacturing Company

36-04 Skillman Avenue

Long Island City, NY 11101

Tel. 718-392-6650