



# Propper Indicator Tapes

## For use in Steam or Ethylene Oxide gas sterilizers

Propper indicator tapes are intended for use in securing packaging materials used in sterilization processes. The printed indicator inks on the tapes are process indicators. They are intended to demonstrate exposure to a sterilization process. They should not be interpreted as indicators of the adequacy of the sterilization cycle.

### **STEAM** Strate-Line<sup>®</sup> Autoclave Tape - Reorder Nos.: 268004, 268005 & 268006.

Universal application: suitable for use with nonwoven fabric wraps & other sterilization packaging materials.

The light colored indicator stripes turn dark when exposed to saturated steam.

**Strate-Line<sup>®</sup> Autoclave Tapes are AAMI/ISO 11140-1, Class 1 Indicators.**

### **STEAM** Blue Autoclave Tape - Reorder Nos.: 268008, 268009.

Developed specifically for greater adhesion. For use with nonwoven disposable wraps only. Do not use on linen or muslin.

The light colored indicator stripes turn dark when exposed to steam sterilization.

**Blue Autoclave Tapes are AAMI/ISO 11140-1, Class 1 Indicators.**

### **EO** Gas-Chex<sup>®</sup> EO Gas Tape- Reorder Nos.: 268015.

Suitable for nonwoven and fabric wraps. The word "gas" changes from red to green when exposed to ethylene oxide.

**NOT MADE WITH  
NATURAL RUBBER LATEX**

Propper Manufacturing Co., Inc.  
36-04 Skillman Ave • Long Island City, NY 11101  
Tel. (718) 392-6650 • [www.proppermfg.com](http://www.proppermfg.com)  
email: [customerservice@proppermfg.com](mailto:customerservice@proppermfg.com)