

Propper Manufacturing Company

Microscope Slides and Cover Glass

Propper Micro Slides and Cover Glass are made of exceptional quality water white glass. They are available in all the most popular sizes for use with blood smears, fluorescence studies and all routine microscopy.

- Clear, non-fogging, non-corrosive first quality water white glass
- Pre-cleaned and ready-to-use
- Designed for uniform wetting
- Guaranteed exceptional clarity and uniformity
- Excellent for long-term storage of samples
- Frosted portion provides an excellent writing surface and opacity



- Premium Micro Slides have precisely cut square shape and unmatched surface quality. The slides are ideal for uniform blood smears and all other routine applications.
- **Slides Size:** 3" x 1" (75.4 x 25.1mm) and 0.95 - 1.05mm thickness
- **Package:** 72 per box, 20 boxes/case
- Bev-L-Edge Micro Slides have smooth ground edges and clipped corners for added safety. The slides are excellent for all manual and automated processes in microscopy.
- **Slides Size:** 3" x 1" (75.5 x 25.5mm) and 0.95 - 1.05mm thickness
- **Package:** 72 per box, 20 boxes/case

Reorder No.	Product
15300100	Bev-L-Edge® Micro Slides Plain Slides
15400100	Bev-L-Edge® Micro Slides Frosted on Two Sides (Twin-Frost)
15810100	Premium Micro Slides Plain Slides
15812100	Premium Micro Slides Frosted on Two Sides (Twin-Frost)



Propper Manufacturing Company, Inc | marketing@proppermfg.com | www.proppermfg.com

36-04 Skillman Ave., Long Island City, NY 11101 | Tel : 718-392-6650 | Fax : 718-482-8909

Propper Manufacturing Company

Microscope Slides and Cover Glass

- Propper Select Cover Glass are for use in all types of microscopy.
- Designed with high resistance to breakage
- Parallel surfaces designed for uniform wetting and homogeneity



Reorder No.	Size	Thickness
14110200	22 x 22mm	0.13-0.17mm
14110700	22 x 40mm	0.13-0.17mm
14111400	24 x 50mm	0.13-0.17mm
14120200	22 x 22mm	0.17-0.25mm
14130200	22 x 22mm	0.16-0.19mm
14130900	22 x 55mm	0.16-0.19mm



Propper Manufacturing Company, Inc | marketing@proppermfg.com | www.proppermfg.com

36-04 Skillman Ave., Long Island City, NY 11101 | Tel : 718-392-6650 | Fax : 718-482-8909