

Propper Laryngoscopes are designed for Longevity, Ease-of-Use and Infection Control!

The Laryngoscope System is

- **Versatile** : 3 handle sizes available with two light sources and 12 blade varieties.
- **Easy To Sterilize** : Handles feature removable inner compartments allowing hassle free cleaning and sterilization. Blades are designed with integrated fiber optic bundles and select quality stainless steel, to eliminate crevices and prevent fluid entrapment.
- **Standardized** : Blades and Handles are Green System compliant.

Stubby Handle



Reorder No.

19958100
Stubby LED

19957800
Stubby Xenon

"C" Handle



Reorder No.

19958200
C LED

19915000
C Xenon

"AA" Handle



Reorder No.

19958300
AA LED

19915200
AA Xenon

Laryngoscope Handles

- Standard Battery Sizes
- Xenon or Superlume LED bulbs
- Ball and lock system and removable inner compartments prevent fluid entrapment allowing hassle free disinfection and sterilization
- Sterilize Steam, EO and Vaporized Hydrogen Peroxide

Superlume® LED

- Brightest output (400% brighter than Xenon-Halogen)
- High color temperature for increased visibility
- LED light source lasts 1,000x longer than a Halogen bulb
- LED's long lifetime will make frequent bulb replacements a thing of the past

Blades, bulbs and handles are guaranteed against manufacturing and material defects for life!

Propper Laryngoscope Blades

Propper Manufacturing blades feature an integrated, high-performance fiber optic bundle. The large diameter fiber optics, combined with state-of-the-art Superlume LED handles, provides an unmatched visualization for laryngoscope users.

Propper Laryngoscope blades are made of select quality 304 Stainless Steel with flawless pinhole free Nickel/Chrome plating. The epoxy sealed integrated fiber optic channel is heat resistant to 392°F. This streamlined system minimizes fluid entrapment and enables sterilization in Steam, EO and Vaporized Hydrogen Peroxide.

Propper blades feature an easy clip system and are Green System compliant.

Propper Laryngoscope Blades' light output will not fall below 2,000 lux for five years when used with Propper LED handles.

Macintosh Blade



Miller Blade



Size	Description	Reorder No.
1	Macintosh Blade	19916100
2	Macintosh Blade	19914400
3	Macintosh Blade	19914500
4	Macintosh Blade	19914600
5	Macintosh Blade	19916200
0	Miller Blade	19915600
1	Miller Blade	19915700
1.5	Miller Blade	19915500
2	Miller Blade	19915800
3	Miller Blade	19915900
4	Miller Blade	19916000

