

Once-A-Day® Bowie & Dick Autoclave Test Record

Once-A-Day® Bowie and Dick Autoclave Test Record



Propper's Once-A-Day® Autoclave Test Record A4 test sheets are designed for use in standard towel packs. Bowie and Dick tests provide a challenge for pre-vacuum sterilizers at the start of each day or after routine or other maintenance.

Once-A-Day® Bowie and Dick Autoclave Test Record features specifically formulated highly reactive ink and easy interpretation of color change guarantees the earliest possible warning of impending problems in pre-vacuum sterilizers.

Once-A-Day® Test Sheet Features:

- Unique printing process gives uniformly reactive ink imprint for decisive color change after testing
- A4 size (210mm x 297mm) size individual Bowie and Dick test sheets meet AAMI ST-79 recommendations and ISO 11140-4 and 11140-5 requirements
- Stable white to black color change after exposure allows for retention of the test sheet for record keeping as long as necessary
- Test sheets detect air leaks and verifies that the sterilizer effectively removes air from the sterilization chamber

Reorder No.

Packaging

26600500

1 Box of 50 Bowie and Dick Test Sheets



Propper Manufacturing Company, Inc | marketing@proppermfg.com | www.proppermfg.com

36-04 Skillman Ave. Long Island City, NY, USA 11101 | Tel : +1-718-392-6650 | Fax : +1-718-482-8909