

O.K.® Sterilizer Bags



Quality and Reliability

Made from 100% medical grade paper, ensuring a sterile barrier

Confidence

O.K.® Bags feature a pre-printed process indicator to allow for easy visualization and confirmation of exposure to steam

Safety

O.K.® Bags are easily secured by folding over the top edge and then sealing with autoclave tape

O.K.® Steam Sterilization Bags

For use in Steam Sterilizers

O.K.® Steam Sterilization Bags are an ideal solution to contain small instruments or implants inside a container when paper-plastic pouches are not appropriate (AAMI ST-79:2017).

All paper design allows for adequate air removal and steam penetration and drying while simultaneously providing a good sterile barrier.

All sizes feature a gusset which allows for bulky objects to be placed safely inside the bag, avoiding breakage.

Features

- Superior quality medical grade paper with reinforced folded bottom maintains the integrity of the bag, even in the case of wet steam or heavy instruments
- The vertical seam is double bonded for added security
- All-paper bag will not become dry or brittle, allowing each item to be stored for maximum recommended period of time
- Surfaces have extra-smooth, machine-glazed finish to eliminate lint
- Process indicator exhibits a pink to black color change on exposure to steam
- Easily seal bags with sterilization tape
- Single use only



Reorder No.	Bags/Box	Length	Width	Height
02100200	1000	2.5in (6.4 cm)	1.5in (3.8cm)	8in (20.3cm)
02100300	1000	2.5in (6.4 cm)	11.5in (3.8cm)	10in (25.4)
02101100	1000	5.75in(13.7cm)	1.5in (3.8cm)	12in (30.5cm)

