


STEAM•DOT™ BLU® PROCESS INDICATORS

Lead & Heavy Metal Free Self-Adhesive Process Indicator for Steam Sterilizers

Environmentally friendly technology and top quality performance

- For use in gravity 250°F and pre-vacuum 270°F, 273°F, and 275°F cycles
- Distinct **blue** to **pink** color change on exposure to steam
- Immediate differentiation between processed and unprocessed packs

Reorder No.	Description	Packaging
36400200	1/2" (1.3cm) Diameter	1000/Roll in dispenser box




STEAM•DOT™ PROCESS INDICATORS

Self-Adhesive Process Indicator for Steam Sterilizers

Fast and positive visual check on items that were sterilized by steam

- For use in gravity 250°F and pre-vacuum 270°F, 273°F, and 275°F cycles
- Distinct **white** to **black** color change on exposure to steam
- Immediate differentiation between processed and unprocessed packs

Reorder No.	Description	Packaging
36400100	3/4" (1.9cm) Diameter	1000/Roll in dispenser box



STERI•DOT® PROCESS INDICATORS

Self-Adhesive Process Indicator for Ethylene Oxide Sterilizers

Fast and positive visual check on items that were sterilized by ethylene oxide

- For use in all ethylene oxide sterilization cycles
- Distinct **red** to **green** color change on exposure to ethylene oxide gas
- Immediate differentiation between processed and unprocessed packs

Reorder No.	Description	Packaging
36100100	9/16" (1.5cm) Diameter	2000/Roll in dispenser box
36100200	1" (2.5cm) Diameter	1000/Roll in dispenser box

