



FDA Approved

Infrared Digital Thermometer

Non-Contact Forehead Measurement
Highly Precise Infrared Sensor
Easy to Clean and Disinfect
Easy-to-Read LED
Digital Display
Temperature Displays
in 1 second

Reorder Number:
26915400



Features

- ◆ **Fever warning signal**
- ◆ **Memory Recall 32 temperatures**
- ◆ **Use 2 AAA batteries**
- ◆ **Auto power-off after 60 seconds**

Product Specifications

Measurement Units : Celsius (°C)/Fahrenheit (°F)
Storage Conditions : -20°C -+55°C (-4°F to 131°F)
with a relative humidity of 15%-93%
Measurement Distance : 15mm to 50mm from the measuring point
Measurement Range : 32.0°C-42.2°C (89.6°F-108.0°F)
Measurement Accuracy : ±0.2°C (35.0°C-42.0°C) \ ±0.4°F (95.0°F-107.6°F);
±0.3°C (±0.5°F) other range
Dimensions : 144mm x 101mm x 36mm
Weight : 88g without batteries
Product Life : 5 years