


Seracult® and Seracult Plus®

Seracult - Designed for use in routine physical examinations in clinical and laboratory settings, Seracult is a qualitative aid in the diagnosis of colorectal cancer.

Seracult Plus - Seracult Plus was developed with your most high risk patients in mind, and can detect smaller concentrations of occult blood in the sample and is a functional equivalent to iFOBT.



Easy and Simple Application



Convenient USPS Approved mailing envelopes



Performance control on every slide

- Simple:** Apply sample to front side of guaiac paper, apply developer to back side of paper, get results in less than 60 seconds if blood is present
- Sensitive:** Seracult detects hemoglobin to 0.38 mg Hb/100ml concentration
Seracult Plus detects hemoglobin to 0.32 mg Hb/100ml and is sensitive to upper GI bleeding
- Secure:** Hydrophobic wax barrier between Performance Control area and sample
- Stable:** 28 day slide sample stability

Seracult and Seracult Plus provide options for all of your needs

- Single Slide:** Optimal for fast results, in all clinical settings, contains performance control and requires sample from one bowel movement
- Triple Slide:** Optimized for at home use, contains performance control, and requires sample from three bowel movements
- Tape:** Best for inpatient use, requires sample from single bowel movement



Seracult® and Seracult Plus®

Reorder No.	Unit	Qty/Unit	Description
37100100	Box	100 Slides	Seracult® Single Slide Fecal Occult Blood Test: Seracult® Single Slides: two 15 ml bottles of developer and 100 applicators
37100200	Case	1000 Slides	Seracult® Single Slide Fecal Occult Blood Test: Seracult® Single Slides: 20 15ml bottles of developer and 1000 applicators
3710025P	Case	1000 slides	Seracult® Single Slide Fecal Occult Blood Test: Seracult® Single Slides
37200300	Case	100 Kits	Seracult® Fecal Occult Blood Test Patient Kit: Seracult Patient Kits: two boxes, each containing 50 kits (triple slide and 3 applicators in dispensing envelope) and three 15 ml bottles of developer
37200400	Box	34 Triple Slides	Seracult® Triple Slide Fecal Occult Blood Test: Seracult Triple Slides: two 15 ml bottles of developer and 102 applicators
37200500	Case	340 Triple Slides	Seracult® Triple Slide Fecal Occult Blood Test: Seracult Triple slides: twenty 15 ml bottles of developer and 1020 applicators
37200700	Case	80 Kits	Seracult® Fecal Occult Blood Test Mailing Kit: 80 mailing kits (triple slide & 3 applicators, 3 collection tissues and USPS approved foil lined mailing envelope) and two 15 ml bottles of developer
37400100	Box	100 Slides	Seracult® Plus Single Slide Fecal Occult Blood Test: two 15ml bottles of developer, 100 applicators
37400200	Case	1000 Slides	Seracult® Plus Single Slide Fecal Occult Blood Test: twenty 15ml bottles of developer, 1000 applicators
37500300	Case	100 Kits	Seracult® Plus Fecal Occult Blood Test Patient Kit: 100 Patient Kits, envelope contains triple slide and three applicators. Each of the two boxes has 50 kits and three 15ml bottles of developer
37500400	Box	34 Triple Slides	Seracult® Plus Triple Slide Fecal Occult Blood Test: 34 Triples Slides, two 15ml bottles of developer, 102 applicators
37500700	Case	80 Kits	Seracult® Plus Fecal Occult Blood Test Mailing Kit: 80 Mailing Kits, envelope contains triple slide, three applicators, three flush-able collection tissues, U.S. Postal Service approved Mailing Envelope (foil lined).
37701000	Box	2 Tape Kits	Seracult® Plus Fecal Occult Blood Test Tape & Developer Kit : 2 rolls of tape 2 bottles developer
37701500	Box	25 Bottles	Seracult® Plus Fecal Occult Blood Test Developer Solution: Twenty 15ml bottles of developer
37901000	Box	2 Tape Kits	Seracult® Fecal Occult Blood Test Tape & Developer Kit : 2 rolls of tape 2 bottles developer
37901600	Box	20 bottles	Seracult® Fecal Occult Blood Test Developer Solution: Seracult Developer: 20 bottles, 15 ml each
37901700	Box	100 Patient Kits	Seracult® Fecal Occult Blood Test Patient Instruction Envelope: 100 Seracult Primary Envelopes with patient instructions
37901500	Box	20 bottles	Seracult® Fecal Occult Blood Test Developer Solution: Seracult® Developer: 20 bottles, 15 ml each

