



Propper Solid Glass Beads

Solid Glass beads are used to reduce bumping during boiling, refluxing, and other distilling processes in laboratories. Propper beads are used in absorption columns and distillers where a large surface of glass is desirable.

- Manufactured from soda lime glass.
- Durable and will not disintegrate or affect delicate compounds.
- Highly resistant to attack by most cold and heat.
- Packaged in approximately 1 Lb zip-lock bag for your convenience.

1mm Diameter

Reorder#: 03000200

Ideal for grinding and dispersing media, paints, pigments and other chemicals and minerals.

May be used as filler material in medical and industrial laboratories.

Reorder No.	Size	Size Range (mm)	Quantity per Lb
03000200	1mm Beads	1.0–1.3	228,000
03000400	2mm Beads	1.8–2.2	43,350
03000600	3mm Beads	2.7–3.3	13,300
03000800	4mm Beads	3.7–4.3	5,420
03000900	5mm Beads	4.7–5.3	2,770
03001000	6mm Beads	5.7–6.3	1,600

2mm ▪ 3mm ▪ 4mm ▪ 5mm ▪ 6mm Diameter

Reorder#: 03000400, 03000600, 03000800, 03000900, 03001000

Best for use for stirring and mixing media in aerosol sprays, as valves for bottle closures, and for the mechanical coating of small metallic pieces, through mechanical plating.

May be used in water wells, to replace filter gravel in the annular space.



Propper Manufacturing Company, Inc | marketing@proppermfg.com | www.proppermfg.com

36-04 Skillman Ave., Long Island City, NY 11101 | Tel : 718-392-6650 | Fax : 718-482-8909