



propper
Manufacturing Co., Inc.

Contact us to order Steam Chex today!
customerservice@proppermfg.com
Tel. (718) 392-6650

Steam Chex is a pressure-sensitive process indicator tape designed to seal packs that are exposed to steam sterilization.

- Made of a saturated crepe paper printed with non-lead based ink, and coated with solvent free adhesive.
- Not made with lead based ink or natural latex rubber.
- Steam Chex features a green to dark brown/black color change.

Reorder No.: Quantity:
26810800 48 Rolls • ¾" x 60 YDS
26810900 36 Rolls • 1" x 60 YDS

