

Steri-Wrap® I

Sterilization Wrap

Product Description

Steri-Wrap I is a blue, water and alcohol repellent wrap. It is made of wet formed nonwoven fabric, manufactured from natural wood pulp and polyester fibers Steri-Wrap I does not contain any natural latex rubber. It can be used to contain instruments during steam, ethylene oxide or gamma radiation sterilization.

Steri-Wrap I should not be used in dry heat, chemical vapor or hydrogen peroxide sterilization.

Steri-Wrap I is available in sheets

31101200 Steri-Wrap® I - 12" x 12" Wrap - 1 Case of 1000
31101500 Steri-Wrap® I - 15" x 15" Wrap - 1 Case of 1000
31102400 Steri-Wrap® I - 24" x 24" Wrap - 1 Case of 500
31103620 Steri-Wrap® I - 36" x 36" Wrap - 1 Case of 200
31104800 Steri-Wrap® I - 48" x 48" Wrap - 1 Case of 100

Benefits

Steri-Wrap I is designed to be low linting and fluid repellent. Steri-Wrap I is best used for containing small or light packages for sterilization. Steri-Wrap I is compliant with CFR 1610 Class I flammability requirements.

Technical Data

Property	Steri-Wrap® I	Unit
Dry Mullen Burst Strength	1725	g/cm ²
Dry Grab Tensile Strength MD	8500	grams
Mason Jar	120	minutes
Alcohol Repellency	7	rating
Air Permeability	80	L/min/100cm ²
Flammability	5.5	seconds
Basis Weight	51-53	gsm
Elmendorf Tear MD	65	g
Suter Waterhead	370	mm



Shelf Life

60 months.

Storage Conditions

Store at temperatures below 32°C (90°F) and relative humidity below 70%.