

Steri-Wrap® II

Sterilization Wrap

Product Description

Steri-Wrap II is a green drapeable all-purpose, water and alcohol repellent fabric wrap. Steri-Wrap II can be used to contain instruments during steam, ethylene oxide or gamma radiation sterilization. It has been manufactured from cellulosic and polyester fibers bonded with an acrylate ester copolymer binder. It does not contain any natural rubber latex.

Steri-Wrap II should not be used in dry heat, chemical vapor or hydrogen peroxide sterilization.

Steri-Wrap II is available in sheets

31201500 Steri-Wrap® II - 15" x 15" Wrap - 1 Case of 1000

31202400 Steri-Wrap® II - 24" x 24" Wrap - 1 Case of 500

31203040 Steri-Wrap® II - 30" x 30" Wrap - 1 Case of 400

31203620 Steri-Wrap® II - 36" x 36" Wrap - 1 Case of 200

31205400 Steri-Wrap® II - 54" x 54" Wrap - 1 Case of 100

Benefits

Steri-Wrap II features excellent bacterial barrier properties, fluid repellency and strength properties. It has been designed as a sterilization packaging material but can also be used as a sterile field. Steri-Wrap II conforms to ISO 11607-1:2019 and EN 868-2:2017.

Technical Data

Property	Steri-Wrap®II	Unit
Dry Mullen Burst Strength	2000	g/cm ²
Dry Grab Tensile Strength MD	8700	grams
Mason Jar	130	minutes
Alcohol Repellency	9	rating
Air Permeability	80	L/min/100cm ²
Flammability	6.0	seconds
Basis Weight	57-58	gsm
Elmendorf Tear MD	125	g
Suter Waterhead	340	mm



Shelf Life

60 months.

Storage Conditions

Store at temperatures below 32°C (90°F) and relative humidity below 70%.



36-04 Skillman Avenue
Long Island City, NY 11101
Tel. (718) 392-6650 • Fax. (718) 482-8909
www.proppermfg.com