

Steri-Wrap® III

Sterilization Wrap

Product Description

Steri-Wrap III is a blue water, alcohol and saline resistant wrap. Steri-Wrap III is made of a wet formed non-woven fabric, manufactured from cellulosic and polyester fibers bonded with an acrylate ester copolymer binder. It does not contain any natural rubber latex. Steri-Wrap III can be used to contain instruments during steam, ethylene oxide or gamma radiation sterilization.

Steri-Wrap III should not be used in chemical vapor, dry heat or hydrogen peroxide sterilization.

Steri-Wrap III is available in sheets

31301800 Steri-Wrap® III - 18" x 18" Wrap - 1 Case of 400
31302440 Steri-Wrap® III - 24" x 24" Wrap - 1 Case of 400
31303600 Steri-Wrap® III - 36" x 36" Wrap - 1 Case of 150
31304810 Steri-Wrap® III - 48" x 48" Wrap - 1 Case of 100
31305400 Steri-Wrap® III - 54" x 54" Wrap - 1 Case of 50

Benefits

Steri-Wrap III is made of low linting fabric and has a good combination of barrier and breathability which makes it highly suited for sterilization and apparel applications. Steri-Wrap III conforms to ISO 11607-1:2019 and EN 868-2:2017.

Technical Data

Property	Steri-Wrap®III	Unit
Dry Mullen Burst Strength	1930	g/cm ²
Dry Grab Tensile Strength MD	9230	grams
Mason Jar	120	minutes
Alcohol Repellency	7	rating
Air Permeability	80	L/min/100cm ²
Flammability	6.4	seconds
Basis Weight	59-61	gsm



Shelf Life

60 months.

Storage Conditions

Store at temperatures below 32°C (90°F) and relative humidity below 70%.



36-04 Skillman Avenue
Long Island City, NY 11101
Tel. (718) 392-6650 • Fax. (718) 482-8909
www.proppermfg.com