

Choose from a manufacturer that offers:
Safety, Versatility and Quality

Steri-Wrap® Single Use, Non-Woven Fabric Sterilization Wrap is available in three grades and variety of sizes. Validated in the most commonly used sterilants making it ideal for use in clinical and industrial settings.

All Steri-Wrap® grades feature:

- **Efficient microbial barrier** to increase pack security
- Tightly anchored fiber structures make wrap **soft, strong and low lint**
- **Excellent liquid repellency** to protect contents from contamination
- **Validated for use in most commonly used sterilization cycles:** pre-vacuum and gravity steam cycles, ethylene oxide, gamma radiation and electron beam radiation.

Three grades of wrap fabric to meet all of your needs

Steri-Wrap® I: Light to medium weight, ideal for general hospital use

Steri-Wrap® II: Medium weight, greater drapeability

Steri-Wrap® III: Heavy weight, suitable for your heaviest instrument packages




Steri-Wrap® I




Steri-Wrap® II



Steri-Wrap® III



Tape will not cause delamination on either side of wrap



Excellent liquid repellency

Size		Steri-Wrap® I		Steri-Wrap® II		Steri-Wrap® III	
Inches	Cm	Reorder No.	Qty/Case	Reorder No.	Qty/Case	Reorder No.	Qty/Case
12x12	31x31	31101200	1000	-	-	-	-
15x15	38x38	31101500	1000	31201500	1000	-	-
18x18	46x46	-	-	-	-	31301800	400
24x24	61x61	31102400	500	31202400	500	31302440	400
30x30	76x76	-	-	31203040	400	-	-
36x36	91x91	31103620	200	31203620	200	31303600	150
45x45	114x114	-	-	-	-	-	-
48x48	122x122	31104800	100	-	-	31304810	100
54x54	137x137	-	-	31205400	1000	31305400	50

