

Propper pressure sensitive sterilization tape is designed to secure sterile packages and indicate exposure to sterilization conditions.

All Propper tape features lead free ink.

### Steam Chex™ Autoclave Indicator Tape

*For use in Steam Sterilizers*

Pressure-sensitive tape featuring strong adhesive for use with disposable wraps during steam sterilization. Steam Chex™ features strong, pressure-sensitive adhesive designed to use on disposable sterilization wraps, like Propper Steri-Wrap® and Steri-Wrap® Unite™. Made without natural latex rubber. Write directly on tape or use Duo-Record® Record Keeping Labels.

**Green to Dark Brown/Black indicator change after exposure**




| Reorder No. | Rolls/Box | Width        | Length/Roll | Core Diameter |
|-------------|-----------|--------------|-------------|---------------|
| 26810800    | 48        | 3/4in (18mm) | 60yd (55m)  | 3in (7.62cm)  |
| 26810900    | 36        | 1in (25mm)   | 60yd (55m)  | 3in (7.62cm)  |

### Verline™ Autoclave Indicator Tape

*For use in Steam Sterilizers*

Verline™ features strong, pressure-sensitive adhesive designed to use on disposable sterilization wraps, like Propper Steri-Wrap® and Steri-Wrap® Unite™. Made without natural latex rubber. Write directly on tape or use Duo-Record® Record Keeping Labels.

**Green to Dark Brown/Black indicator change after exposure**




| Reorder No. | Rolls/Box | Width        | Length/Roll | Core Diameter |
|-------------|-----------|--------------|-------------|---------------|
| 26810500    | 48        | 3/4in (18mm) | 60yd (55m)  | 3in (7.62cm)  |
| 26810600    | 36        | 1in (25mm)   | 60yd (55m)  | 3in (7.62cm)  |

### Gas-Chex® EO Gas Indicator Tape

*For use with Ethylene Oxide Gas Sterilizers*

Universal adhesive for all types of reusable and disposable wrap, paper, paper/plastic blends metal and glass. Write directly on tape or use Duo-Record® Record Keeping Labels.

**Word "gas" exhibit a Red to Green color change on exposure to EO**



| Reorder No. | Rolls/Box | Width        | Length/Roll | Core Diameter |
|-------------|-----------|--------------|-------------|---------------|
| 26801500    | 48        | 3/4in (18mm) | 60yd (55m)  | 3in (7.62cm)  |

