



**Strate-Line® is no longer available for purchase.
Contact us to order Verline™ today!
customerservice@proppermfg.com
Tel. (718) 392-6650**



Strate-Line and Verline Autoclave Indicator Tapes are pressure-sensitive process indicator tapes designed to seal packs that are exposed to steam sterilization.



Strate-Line and Verline Autoclave Indicator Tapes are made of saturated crepe paper and coated with solvent-free adhesive.



Strate-Line and Verline Autoclave Indicator Tapes are not made with lead based ink or natural latex rubber.



Strate-Line Autoclave Indicator Tape features a WHITE to dark brown/black color change.

Verline Autoclave Indicator Tape features a GREEN to dark brown/black color change.

Reorder No.:


26810500

26810600

Quantity:

48 Rolls • ¾" x 60 YDS

36 Rolls • 1" x 60 YDS



@proppermfg