



Emu-Chex™ Instructions for Use

Instructions for Use for 121°C 20 Minute Emu-Chex™ Indicator

Emu-Chex™ is emulating indicator designed to react to all the critical variables of the steam sterilization process. The 121°C, 20 min indicator can be used in either gravity or pre-vacuum cycles.

Place an emulating chemical indicator into each pack, pouch or tray for sterilization in the most difficult area for steam to reach. For large packs, place Emu-Chex™ at several locations. Run cycle. After processing, when sterile presentation is made, examine the indicator dot. A brown color indicates a pass result. Grey, blue, green, black and mottled colors indicator underprocessing: do not use the contents of the pack. Determine and correct the cause of failure and reprocess the content according to your hospital procedures.

Instructions for Use for 132°C 4 Minute Indicator

Emu-Chex™ is emulating indicator designed to react to all the critical variables of the steam sterilization process. The indicator can be used in pre-vacuum cycles.

Place an emulating chemical indicator into each pack, pouch or tray for sterilization in the most difficult area for steam to reach. For large packs, place Emu-Chex™ at several locations. Run cycle. After processing, when sterile presentation is made, examine the indicator dot. A brown color indicates a pass result. Grey, blue, green, black and mottled colors indicator underprocessing: do not use the contents of the pack. Determine and correct the cause of failure and reprocess the content according to your hospital procedures.

Instructions for Use for 135°C 3 Minute Indicator

Emu-Chex™ is emulating indicator designed to react to all the critical variables of the steam sterilization process. The 1 indicator can be used in either gravity or pre-vacuum cycles.

Place an emulating chemical indicator into each pack, pouch or tray for sterilization in the most difficult area for steam to reach. For large packs, place Emu-Chex™ at several locations. Run cycle. After processing, when sterile presentation is made, examine the indicator dot. A brown color indicates a pass result. Grey, blue, green, black and mottled colors indicator underprocessing: do not use the contents of the pack. Determine and correct the cause of failure and reprocess the content according to your hospital procedures.