


Solid Soda Lime Glass Beads



Propper Solid Chemical Glass Beads are designed for use as mixing beads, to reduce bumping during boiling, refluxing and other distilling processes or for packing distillation columns.

- Soda Lime Glass Beads are durable and will not disintegrate
- Propper glass beads are highly resistant to attack by most cold or hot liquids and vapors
- For use wherever a large surface area of glass is desirable

Packaged in 1 pound boxes for your convenience.

Reorder No.	Size (mm)	Approx No/lb
030002	1.0	400,000
030004	2.0	50,000
030006	3.0	12,000
030008	4.0	4,500
030009	5.0	2,500
030010	6.0	1,500

