

Microscope Slides and Cover Glass

Propper Microscope Slides and Cover Glass are made of exceptional quality water-white glass. They are available in all the most popular sizes for use with all levels of microscopy.

- Clear, non-fogging, non-corrosive first quality water-white glass
- Ready-to-use: Cleaned and Washed
- Designed for uniform wetting
- Guaranteed exceptional clarity and uniformity
- Precision ground edges



- Propper Premium Micro Slides are available in plain or frosted styles
- Frosted Style provides an excellent writing surface and opacity
- Slides Size: 3"x1" (76.2x25.4mm) and 0.9-1mm thickness

Reorder No.	Plain	Package
158101		72/box, 20/case
	Frosted on Two Sides:Twin-Frost	
158121		72/box, 20/case



- Bev-L-Edge Micro Slides have smooth ground bevelled edges for uniform smears and are perfect for blood smears and routine microscopy
- Available in plain and frosted styles
- Slides Size: 3"x1" (76.2x25.4mm) and 0.9-1mm thickness

Reorder No.	Plain	Package
153001		72/box 20/case
	Frosted on Two Sides:Twin-Frost	
154001		72/box 20/case



- Propper Select Cover Glass are for use in all types of microscopy.
- Designed with high resistance to breakage
- Parallel surfaces designed for uniform wetting and homogeneity

Reorder No.	Size	Thickness
141102	22x22	0.13-0.17mm
141302	22x22	0.16-0.19mm
141202	22x22	0.17-0.25mm
141107	22x40	0.13-0.17mm
141309	22x50	0.16-0.19mm
141114	24x50	0.13-0.17mm

