

# Sterilization Packaging: Dust Covers

## Propper Dust Covers


Propper self-seal dust covers protect sterile packs from dust, moisture, physical damage, or airborne contaminants during storage, handling, and transportation.

Both self-seal and nonslip dust covers are made of 76 $\mu$  plastic providing superior strength and quality. All Propper Dust Covers comply with AAMI ST-79:2017

### Self-Seal Dust Covers Features:

- **Self sealing:** Eliminates the need for heat sealers
- **Wide adhesive band:** Provides an easy sealing procedure
- **Quality:** 76 $\mu$  plastic provides excellent strength
- **Easy opening:** Linear tear plastic with notches on both sides
- **Easy labeling:** Special white area for recording contents and date
- **Clarity:** Transparent on both sides for easy content identification

Reorder No.	Bags/Case	Inches	Centimeters
32121200	1000	8x12	20.3x30.5
32121500	1000	10x15	25.4x38.1
32122000	500	12x20	30.5x50.8
32122200	500	16x22	40.6x55.9
32123000	250	16x30	40.6x76.2
32123200	250	22x30	55.9x76.2



### Nonslip Dust Covers Features:

- **Ease:** Scored line facilitates easy opening
- **Strength:** 76 $\mu$  plastic exceeds the industry standard
- **Nonslip plastic:** Prevents sterilized packs from sliding off delivery to operating room or procedure room
- **Record keeping:** Specially treated matte area for writing date, number, content info and operators initials
- **Secure:** Covers are easily sealed with tape

Reorder No.	Bags/Case	Inches	Centimeters
32101500	1000	10x15	25.4x38.1
32102000	500	12x20	30.5x50.8
32102200	500	16x22	40.6x55.9
32103000	250	16x30	40.6x76.2
32103200	250	22x30	55.9x76.2

