

Propper Manufacturing Company

Steam & Ethylene Oxide Process Indicators

Steam Dot® Blu™

Lead & Heavy Metal Free Self Adhesive Process Indicator for Steam Sterilizers



Environmentally Friendly Technology and Top Quality Performance

- For use in cycles at 250F, 270f, 272F & 275F
- Distinct **blue** to **pink** color change on exposure to steam
- Immediate differentiation between processed and unprocessed packs

Reorder no.	Description	Packaging
36400200	1/2" Steam-Dot® Blu™	1000/roll in dispenser box

Steam Dot® Indicators

Self-Adhesive Process Indicator for Steam Sterilizers

Fast and positive visual check on items that were sterilized by steam

- For use in cycles at 250F, 270f, 272F & 275F
- Distinct **white** to **black** color change on exposure to steam
- Immediate differentiation between processed and unprocessed packs



Reorder no.	Description	Packaging
36400100	3/4" Steam-Dot®	1000/roll in dispenser box

Steri-Dot® Indicators

Self-Adhesive Process Indicator for Ethylene Oxide Sterilizers



Fast and positive visual check on items that were sterilized by Ethylene Oxide

- For use in all Ethylene Oxide Sterilization cycles
- Distinct **red** to **green** color change on exposure to steam
- Immediate differentiation between processed and unprocessed packs

Reorder no.	Description	Packaging
36100100	9/16' Steri-Dot®	2000/roll in dispenser box
36100200	1" Steri-Dot®	1000/roll in dispenser box

