

Cat. Number: 266115
24 packs per case

ONCE·A·DAY[®] ecos[™]

Lead & Heavy Metal Free Bowie & Dick Test for Prevacuum Sterilizers

Eco Friendly

- Non-toxic, no environmental pollutants,
- Contains no lead or other heavy metals
- Latex Free
- 100% Biodegradable
- Optimized test pack size reduces waste

Best-in-Class Performance

- Reveals micro-leaks, residual air & non-condensable gases
- Detects low temperature/pressure conditions, wet & superheated steam
- Guarantees distinct color change
- Processed test sheets will not revert

Full compliance with all requirements of:

ISO 11140-1:2006, Class 2
ISO 11140-4:2007

STEAM (134°-137°C) up to
3.5 minute exposure

36-04 Skillman Ave • Long Island City, NY 11101 • USA
Tel. (718) 392-6650 • Fax: (718) 482-8909 • www.proppermfg.com
European contact: Proplink Medical, Ltd. • Telford, TF2 9HW, UK
tel: +44(0) 1952 610962 • info@proplinkmedical.co.uk

propper
Manufacturing Co., Inc.

