


Propper pressure sensitive sterilization tape is designed to secure sterile packages and indicate exposure to sterilization conditions.

### Steam Chex™ Autoclave Indicator Tape

*For use in Steam Sterilizers*

Pressure-sensitive tape featuring strong adhesive for use with disposable wraps during steam sterilization. Steam Chex™ features strong, pressure-sensitive adhesive designed to use on disposable sterilization wraps, like Steri-Wrap®. Made without natural latex rubber. Write directly on tape or use Duo-Record® Record Keeping Labels.

*Green to Dark Brown/Black indicator change after exposure*




Reorder No.	Rolls/Box	Width	Length/Roll	Core Diameter
26810800	48	3/4in (18mm)	60yd (55m)	3in (7.62cm)
26810900	36	1in (25mm)	60yd (55m)	3in (7.62cm)

### Verline™ Autoclave Indicator Tape

*For use in Steam Sterilizers*

Verline™ features strong, pressure-sensitive adhesive designed to use on disposable sterilization wraps, like Steri-Wrap®. Made without natural latex rubber. Write directly on tape or use Duo-Record® Record Keeping Labels.

*Green to Dark Brown/Black indicator change after exposure*



Reorder No.	Rolls/Box	Width	Length/Roll	Core Diameter
26810500	48	3/4in (18mm)	60yd (55m)	3in (7.62cm)
26810600	36	1in (25mm)	60yd (55m)	3in (7.62cm)

### EO Chex™ Indicator Tape

*For use with Ethylene Oxide Gas Sterilizers*

Universal adhesive for all types of reusable and disposable wrap, paper, plastic, metal and glass. Write directly on tape or use Duo-Record® Record Keeping Labels.

*Pink to Orange/Brown indicator change after exposure*



Reorder No.	Rolls/Box	Width	Length/Roll	Core Diameter
26811500	48	3/4in (18mm)	55yd (50.3m)	3in (7.62cm)
26811600	36	1in (25mm)	55yd (50.3m)	3in (7.62cm)

