

**Clearer Vision, Safer Care:
The Superior Single-Use Choice.**



Propper Manufacturing Company's Magnalume™ plastic disposable laryngoscope blades utilize integrated design with high-performance plastic light transmitter. Combined with state-of-the-art Superlume® LED handles, Propper disposable blades provide an unmatched visualization of the larynx and vocal cords. Magnalume™ blades feature an easy clip system and are ISO 7376 Green System compliant. Magnalume™ disposable blades are packed in a clean environment and sold non-sterile¹.

Intended Use:

Propper disposable laryngoscope blades are medical devices used by healthcare professionals to examine and visualize the larynx and facilitate tracheal intubation.

Key features of Magnalume™ disposable laryngoscope blades:

Easy to use: Disposable blades are ready for use, and don't require cleaning and disinfection.

Better infection control: Disposable blades reduce risk of infection.

Superior strength: Blades are made of superior quality, high strength white Polycarbonate plastic and reinforced arching.

Bright illumination: Light transmission piece is made of clear non-fogging plastic.

Non-magnetic metal components²: Lock system balls, spring, and ball holders are made of low- and non-magnetic metals.

Standardized: Blades are ISO 7376 Green System compliant and fit on any handle compatible with the Green ISO 7376 standard. The green mark is an easily identifiable dot on the right side of each blade.



Propper Manufacturing Company

Magnalume™ Disposable Laryngoscopes Blades

Packaging and Labeling:

- Individual blades are packed in sealed plastic sleeves with labels.
- 10 blades of the same type and size are packed, along with instructions for use, in a cardboard box, which is the unit of sale.

Storage Conditions:

Disposable plastic laryngoscope blades should be stored in a clean, dry, and well-ventilated environment at room temperature in their original packaging. The product must be protected from excessive heat, direct sunlight, humidity, and mechanical damage. Discard after expiration date.

Disposal:

Do not sterilize, reprocess, or re-use. Blades are contaminated after intubation, and it is the user's responsibility to ensure that all used blades are disposed in accordance with medical waste disposal guidelines.

Macintosh Blade



Miller Blade



Size	Description	Reorder No.
1	Macintosh Disposable Laryngoscope Blade	19924100
2	Macintosh Disposable Laryngoscope Blade	19924200
3	Macintosh Disposable Laryngoscope Blade	19924300
4	Macintosh Disposable Laryngoscope Blade	19924400
0	Miller Disposable Laryngoscope Blade	19934000
1	Miller Disposable Laryngoscope Blade	19934100
2	Miller Disposable Laryngoscope Blade	19934200
3	Miller Disposable Laryngoscope Blade	19934300
4	Miller Disposable Laryngoscope Blade	19934400

Note: 1. Magnalume™ blades do not require additional cleaning or disinfection.
2. Blades are not certified for MRI compatibility.

